



Recombinant Feline sarcoma virus Gag polyprotein (gag)

Product Code	CSB-MP366135FAZ
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	P04322
Product Type	Recombinant Protein
Immunogen Species	Feline sarcoma virus (strain Hardy-Zuckerman 4)
Purity	>85% (SDS-PAGE)
Sequence	MSGASSGTAI GAHLFGVSPE CRVLIGDEGA GPSKSLSEVS FSVWYQSRAA RLVIFCLVAS FLVPCLTFLI AETVMGQTIA TPLSLTLDHW SEVRARAHNQ GVEVRKKKWV TLCEAEWVMM NVGWPREGTF SLDNISQVEK KIFAPGPYGH PDQVPYITTW RSLATDPPSW VRPFLPPKP PTPLPQPLSP QPSAPLTSSL YPVLPKTDPP KPPVLPPDPS SPLIDLLTEE PPPYPGGHGP LPSGPRTPTA SPIASRLRER RENPAEESQA LPLREGPNR PQYWPFASD LYNWKSHNPP FSQDPVALTN LIESILVTHQ PTWDDCQQLL QALLTGEERQ RVLLEARKQV PGEDGRPTQL PNVIDETFPL TRPNWDFATP AGREHLRLYR QLLLAGLRGA ARRPTNLAQV KQVV
Source	Mammalian cell
Target Names	gag
Protein Names	Recommended name: Gag polyprotein Alternative name(s): Core polyprotein Cleaved into the following 3 chains: 1. Matrix protein p15 Short name= 2. MA 3. RNA-binding phosphoprotein p12 Alternative name(s): pp12 Capsid protein
Expression Region	1-414
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.