



Recombinant Viral myc transforming protein (MYC)

Product Code	CSB-BP366221ASZ
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	P0C0N8
Product Type	Recombinant Protein
Immunogen Species	Avian retrovirus MH2E21
Purity	>85% (SDS-PAGE)
Sequence	MPLSVSLPSK NYDYDYDSVQ PYFYFEEEEEE NFYLAAQQRS SELQPPAPSE DIWKKFELLP APPLSPSCRS NLAAASCFPS TADQLEMVTE LLGGDMVNQS SICDPDES F VKSIIIRDCM WSGFSAAAKL EKVVSEK LAT YKASRREGGP AAASRPGPPP SGPPPPAGP AASAGLYLHD LGAAAAGCIG SSVVFP CPLG RRGPPGAGPA ALLGVDAPPT AGGGSEEEQE EDEEIDVVT L AEANESESST ESSTEASEEH CKPHHSPLVL ERCHVNIHQH NYAAPPSTKV EYPAKR LKL DSGRVLKQVS NNRKCSPRT SDSEVNDKRR THNVLERQRR NELKLSFFAL RDQIPEVANN EKAPKVVILK RATEYVLSIQ SDEHRLIAEK EQLRRRREQL KHKLEQLRNS RA
Source	Baculovirus
Target Names	MYC
Protein Names	Recommended name: Viral myc transforming protein Short name= v-Myc
Expression Region	1-412
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.