



Recombinant Periplasmic serine endoprotease DegP (degP)

Product Code	CSB-BP366225EOD
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	P0C0V1
Product Type	Recombinant Protein
Immunogen Species	Escherichia coli O157:H7
Purity	>85% (SDS-PAGE)
Sequence	AETS SATTAQQMPS LAPMLEKVMP SVVSINVEGS TTVNTPRMPR NFQQFFGDDS PFCQEGSPFQ SSPFCQGGQG GNGGGQQQKF MALGSGVIID ADKGYVVTNN HVVDNATVIK VQLSDGRKFD AKMVGKDPRS DIALIQQNP KNLTAIKMAD SDALRVGDYT VAIGNPFGLG ETVTSGIVSA LGRSGLNAEN YENFIQTDAA INRGNSGGAL VNLNGELIGI NTAILAPDGG NIGIGFAIPS NMVKNLTSQM VEYGQVKRGE LGIMGTELNS ELAKAMKVDA QRGAFVSQVL PNSSAAKAGI KAGDVITSLN GKPISSFAAL RAQVGTMPVG SKLTLGLLRD GKQVNVNLEL QQSSQNQVDS SSIFNGIEGA EMSNKGKDQG VVVNNVKTGT PAAQIGLKKG DVIIGANQQA VKNIAELRKV LDSKPSVLAL NIQRGDSTIY LLMQ
Source	Baculovirus
Target Names	degP
Protein Names	Recommended name: Periplasmic serine endoprotease DegP EC= 3.4.21.107 Alternative name(s): Heat shock protein DegP Protease Do
Expression Region	27-474
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Full Length of Mature Protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.