



Recombinant Human T-cell leukemia/lymphoma protein 1B (TCL1B)

Product Code	CSB-MP023314HU
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	O95988
Product Type	Recombinant Protein
Immunogen Species	Homo sapiens (Human)
Purity	>85% (SDS-PAGE)
Sequence	MASEASVRLG VPPGRLWIQR PGIYEDEEGR TWVTVVVRFN PSRREWARRAS QGSRYEPSIT VHLWQMAVHT RELLSSGQMP FSQLPAVWQL YPGRKYRAAD SSFWEIADHG QIDSMEQLVL TYQPERKD
Source	Mammalian cell
Target Names	TCL1B
Protein Names	Recommended name: T-cell leukemia/lymphoma protein 1B Alternative name(s): Oncogene TCL-1B Short name= Oncogene TCL1B SYN-1 Syncytiotrophoblast-specific protein TCL1/MTCP1-like protein 1
Expression Region	1-128
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.