



# Recombinant Human Zinc finger BED domain-containing protein 1 (ZBED1)

<b>Product Code</b>	CSB-BP026306HU
<b>Storage</b>	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
<b>Uniprot No.</b>	O96006
<b>Product Type</b>	Recombinant Protein
<b>Immunogen Species</b>	Homo sapiens (Human)
<b>Purity</b>	>85% (SDS-PAGE)
<b>Sequence</b>	MENKSLESSQ TDLKLVAHPR AKSKVWKYFG FDTNAEGCIL QWKKIYCRIC MAQIAYSGNT SNLSYHLEKN HPEEFCEFEVK SNTEQMREAF ATAFSKLKPE SSQQPGQDAL AVKAGHGYDS KKQQELTAAV LGLICEGLYP ASIVDEPTFK VLLKTADPRY ELPSRKYIST KAPEKYGAV REVILKELAE ATWCGISTDM WRSENQNRAY VTAAHFLGL GAPNCLSMGS RCLKTFEVPE ENTAETITRV LYEVFIEWGI SAKVFGATTN YGKDIVKACS LLDVAVHMPC LGHTFNAGIQ QAFQLPKLGA LLSRCRKLVE YFQQSAVAMY MLYEKQKQQN VAHCMLVSNR VSWWGSTLAM LQRLKEQQFV IAGVLVEDSN NHHLMLEASE WATIEGLVEL LQPFKQVAEM LSASRYPTIS MVKPLLHMLL NTTLNKETD SKELSMKEV IAKELSKTYQ ETPEIDMFLN VATFLDPRYK RLPFLSAFER QQVENRVVEE AKGLLDKVVD GGYRPAEDKI FVPPEPPVK KLMRTSTPPP ASVINNMLAE IFCQTGGVED QEEWHAQVVE ELSNFKSQKV LGLNEDPLKW WSDRLALFPL LPKVLQKYWC VTATRVAPER LFGSAANVVS AKRNRLAPAH VDEQVFLYEN ARSGAEAEPE DQDEGEWGLD QEQVFSLGDG VSGGFFGIRD SSFL
<b>Source</b>	Baculovirus
<b>Target Names</b>	ZBED1
<b>Protein Names</b>	Recommended name: Zinc finger BED domain-containing protein 1 Alternative name(s): Putative Ac-like transposable element dREF homolog
<b>Expression Region</b>	1-694
<b>Notes</b>	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
<b>Tag Info</b>	Tag type will be determined during the manufacturing process.
<b>Protein Length</b>	Full length protein
<b>Reconstitution</b>	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
<b>Shelf Life</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life



of lyophilized form is 12 months at -20°C/-80°C.