



Recombinant Human Zinc finger BED domain-containing protein 1 (ZBED1)

Product Code	CSB-MP026306HU
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	O96006
Product Type	Recombinant Protein
Immunogen Species	Homo sapiens (Human)
Purity	>85% (SDS-PAGE)
Sequence	MENKSLESSQ TDLKLVHPR AKSKVWKYFG FDTNAEGCIL QWKKIYCRIC MAQIAYSGNT SNLSYHLEKN HPEEFCEVVK SNTEQMREAF ATAFSKLKPE SSQQPGQDAL AVKAGHGYDS KKQQELTAAV LGLICEGLYP ASIVDEPTFK VLLKTADPRY ELPSRKYIST KAPEKYGAV REVILKELAE ATWCGISTDM WRSENQNRAY VTLAAHFLGL GAPNCLSMGS RCLKTFEVPE ENTAETITRV LYEVFIEWGI SAKVFGATTN YGKDIVKACS LLDVAVHMPC LGHTFNAGIQ QAFQLPKLGA LLSRCRKLVE YFQQSAVAMY MLYEKQKQQN VAHCMLVSNR VSWWGSTLAM LQRLKEQQFV IAGVLVEDSN NHHLMLEASE WATIEGLVEL LQPFKQVAEM LSASRYPTIS MVKPLLHMLL NTTLNKETD SKELSMKEV IAKELSKTYQ ETPEIDMFLN VATFLDPRYK RLPFLSAFER QQVENRVVEE AKGLLDKVVD GGYRPAEDKI FPVPEEPPVK KLMRTSTPPP ASVINNMLAE IFCQTGGVED QEEWHAQVVE ELSNFKSQKV LGLNEDPLKW WSDRLALFPL LPKVLQKYWC VTATRVAPER LFGSAANVVS AKRNRLAPAH VDEQVFLYEN ARSGAEAEPE DQDEGEWGLD QEQVFSLGDG VSGGFFGIRD SSFL
Source	Mammalian cell
Target Names	ZBED1
Protein Names	Recommended name: Zinc finger BED domain-containing protein 1 Alternative name(s): Putative Ac-like transposable element dREF homolog
Expression Region	1-694
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life



of lyophilized form is 12 months at -20°C/-80°C.