



# Recombinant Human Trehalase (TREH)

<b>Product Code</b>	CSB-MP024402HU
<b>Storage</b>	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
<b>Uniprot No.</b>	O43280
<b>Product Type</b>	Recombinant Protein
<b>Immunogen Species</b>	Homo sapiens (Human)
<b>Purity</b>	≥85% (SDS-PAGE)
<b>Sequence</b>	LPPPCES EIYCHGELLN QVQMAKLYQD DKQFVDMPLS IAPEQVLQTF TELSRDHNHS IPREQLQAFV HEHFQAKGQE LQPWTPADWK DSPQFLQKIS DAKLRAWAGQ LHQLWKKLGK KMKPEVLSHP ERFSLIYSEH PFIVPGGRFV EFYYWDSYWV MEGLLLSEMA ETVKGMLQNF LDLVKTYGHV PNGGRVYYLQ RSQPPLLTLM MDCYLHTHTND TAFLQENIET LALELDFWTK NRTVSVSLEG KNYLLNRYVYV PYGGPRPESY SKDVELADTL PEGDREALWA ELKAGAESGW DFSSRWLIGG PNPNSLSGIR TSKLVPVDLN AFLCQAEELM SNFYSLRGND SQATKYRILR SQRLAALNTV LWDEQTGAWF DYDLEKKKKKN REFYPSNLTP LWAGCFSDPG VADKALKYLE DNRILTYQYG IPTSLQKTGQ QWDFPNAWAP LQDLVIRGLA KAPLRRRAQEV AFQLAQNWIR TNFDVYSQKS AMYEKYDVSN GGQPGGGGEY EVQEGFGWTN GVVLMMLLDY GDRLTS
<b>Source</b>	Mammalian cell
<b>Target Names</b>	TREH
<b>Protein Names</b>	Recommended name: Trehalase EC= 3.2.1.28 Alternative name(s): Alpha, alpha-trehalase Alpha, alpha-trehalose glucohydrolase
<b>Expression Region</b>	24-556
<b>Notes</b>	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
<b>Tag Info</b>	Tag type will be determined during the manufacturing process.
<b>Protein Length</b>	Full Length of Mature Protein
<b>Target Details</b>	This gene encodes an enzyme that hydrolyses trehalose, a disaccharide formed from two glucose molecules found mainly in fungi, plants, and insects. A partial duplication of this gene is located adjacent to this locus on chromosome 11.
<b>Reconstitution</b>	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
<b>Shelf Life</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.