



# Recombinant Pig Cholesterol 7-alpha-monooxygenase (CYP7A1), partial

<b>Product Code</b>	CSB-BP006460PI
<b>Storage</b>	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
<b>Uniprot No.</b>	O46491
<b>Product Type</b>	Recombinant Protein
<b>Immunogen Species</b>	Sus scrofa (Pig)
<b>Purity</b>	≥85% (SDS-PAGE)
<b>Sequence</b>	MMSISLLGGI VTAVCCCLWL LLGMRRRQTG EPPLengiIP YLGCALQFGA NPLEFLRANQ RKHGHIFTCQ LMGNYVHFIT NPLSYHKVLC HGKYLDWKKF HFTASAKAFG HRSIDPSDGN TTDNINKTII KTLQGDALNL LAAAMMENLQ LVLRPQVAPQ PEKPAWVTEG MYSFYRVMF EAGYVTLFGK DFIGHDAQKA LILNNLDNFK QFDKIFPALV AGFPIHVFKT GHYAREKLAE GLRLQKLRKR DHISELVRFL NDTLSTLDDA EKAKSLLAVL WASQANTIPA TFWCLFQIR SPEAMKAASE EVNKTLEKAG QKISLDDKPI YLNQIELDSM PVLDSIIKES LRLSSASLNI RTAKEDFTLH LQDGSYNIRK DDIALYPQL MHLDPFIYPD PLTFKYDRYL DENGKTKTTF YSHGLKLYYY YMPFGSGATI CPGRLFVQVE IKQFLILMLS YFDLELVESH VKCPPLDQSR AGLGILPPSN DIEFRYKLRK L
<b>Source</b>	Baculovirus
<b>Target Names</b>	CYP7A1
<b>Protein Names</b>	Recommended name: Cholesterol 7-alpha-monooxygenase EC= 1.14.13.17 Alternative name(s): CYPVII Cholesterol 7-alpha-hydroxylase Cytochrome P450 7A1
<b>Expression Region</b>	1-501
<b>Notes</b>	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
<b>Tag Info</b>	Tag type will be determined during the manufacturing process.
<b>Protein Length</b>	Partial
<b>Target Details</b>	This gene encodes a member of the cytochrome P450 superfamily of enzymes. The cytochrome P450 proteins are monooxygenases which catalyze many reactions involved in drug metabolism and synthesis of cholesterol, steroids and other lipids. This endoplasmic reticulum membrane protein catalyzes the first reaction in the cholesterol catabolic pathway in the liver, which converts cholesterol to bile acids. This reaction is the rate limiting step and the major site of regulation of bile acid synthesis, which is the primary mechanism for the removal of cholesterol from the body.
<b>Reconstitution</b>	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final



concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.

---

**Shelf Life**

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.

Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.