



# Recombinant Mouse NF-kappa-B inhibitor epsilon (Nfkbie)

|                          |   |
|--------------------------|---|
| <b>Product Code</b>      | CSB-EP015764MO-B  |
| <b>Storage</b>           | Store at -20°C, for extended storage, conserve at -20°C or -80°C.   |
| <b>Uniprot No.</b>       | O54910  |
| <b>Product Type</b>      | Recombinant Protein   |
| <b>Immunogen Species</b> | Mus musculus (Mouse)  |
| <b>Purity</b>            | >85% (SDS-PAGE)   |
| <b>Sequence</b>          | MSDARKGPDE ADDSQCDSDGI ESLRSLRSLP EPTAAPGSGS<br>SQSGCPQPWR HAPETHKEPE KEDADGERAD STYASSSLTE SFLLERPEA<br>KDPSPAPGS PLPPAGVLSP QQLEALTYIS EDGDTLLHLA VIHEAPSVLF<br>CCLAFLPQEV LDIQNNLYQT ALHLAVHLDQ PDVVRALVLK GASRILQDQH<br>GDTALHVACR RQNLACACCL LEEQPEPGRQ LSHPLDLQLK NWQGLACLHI<br>ATLQRNQPLI ELLQNGADI DVQEGTSGKT ALHLAVETQE RSLVQFLLQA<br>GARVDARMLN GCTPLHLAAG RGLNSISSTL CEAGADSLLL NVEDETPQDL<br>AEDLLSYLPF DDLKISGKPL LCTD |
| <b>Source</b>            | E.coli  |
| <b>Target Names</b>      | Nfkbie  |
| <b>Protein Names</b>     | Recommended name: NF-kappa-B inhibitor epsilon Short name= NF-kappa-BIE<br>Alternative name(s): I-kappa-B-epsilon Short name= IκB-E Short name= IκB-epsilon Short name= IκappaBepsilon  |
| <b>Expression Region</b> | 1-364   |
| <b>Notes</b>             | Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.   |
| <b>Tag Info</b>          | Tag type will be determined during the manufacturing process.   |
| <b>Protein Length</b>    | Full length protein   |
| <b>Reconstitution</b>    | We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.                         |
| <b>Shelf Life</b>        | The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.<br>Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.  |