



Recombinant Human Vinexin (SORBS3)

Product Code	CSB-BP022406HU
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	O60504
Product Type	Recombinant Protein
Immunogen Species	Homo sapiens (Human)
Purity	≥85% (SDS-PAGE)
Sequence	<p>MQGPPRSLRA GLSLDDFIPG HLQSHIGSSS RGTRVPVIRN GGSNTLNFQF HDPAPRTVCN GGYTPRRDAS QHPDPAWYQT WPGPGSKPSA STKIPASQHT QNWSATWTKD SKRRDKRWVK YEGIGPVDES GMPIAPRSSV DRPRDWYRRM FQQIHRKMPD LQLDWTFEFP PRDPRHLGAQ QRPAHRPGPA TSSSGRSWDH SEELPRSTFN YRPGAFSTVL QPSNQVLRRR EKVDNVWTEE SWNQFLQELE TGQRPKKPLV DDPGEKPSQP IEVLLERELA ELSAELDKDL RAIETRLPSP KSSPAPRRAP EQRPPAGPAS AWSSSYPHAP YLGSARSLSP HKMADGGSPF LGRRDFVYPS STRDPSASNG GGSPARREEK KRKAARLKFD FQAQSPKELT LQKGDIVYIH KEVDKNWLEG EHHGRLGIFP ANYVEVLPAD EIPKPIKPPT YQVLEYGEAV AQYTFKGDLE VELSFRRKGEH ICLIRKVNEN WYEGRITGTG RQGIFPASYV QVSREPLRL CDDGPQLPTS PRLTAAARSA RHPSSPSALR SPADPIDLGG QTSPRRTGFS FPTQEPRPQT QNLGTPGPAL SHSRGPSHPL DLGTSSPNTS QIHWTPYRAM YQYRPQNEDE LELREGDRVD VMQQCDDGWF VGVSRRTQKF GTFPGNYVAP V</p>
Source	Baculovirus
Target Names	SORBS3
Protein Names	Recommended name: Vinexin Alternative name(s): SH3-containing adapter molecule 1 Short name= SCAM-1 Sorbin and SH3 domain-containing protein 3
Expression Region	1-671
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.