



Recombinant Human Rab11 family-interacting protein 3 (RAB11FIP3)

Product Code	CSB-BP019157HU
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	O75154
Product Type	Recombinant Protein
Immunogen Species	Homo sapiens (Human)
Purity	≥85% (SDS-PAGE)
Sequence	MASAPPASPP GSEPPGPDPE PGGPDGPGAA QLAPGPAELR LGAPVGGPDP QSPGLDEPAP GAAADGGARW SAGPAGLEG GPRDPGPSAP PPRSGPRGQL ASPDAPGPGP RSEAPLPELD PLFSWTEEPE ECGPASCPES APFRLQGSSS SHRARGEVDV FSPFPAPTAG ELALEQGP GS PPQPSDLSQT HPLPSEPVGS QEDGPRLRAV FDALDGDGDG FVRIEDFIQF ATVYGAEQVK DLTKYLDPSG LGVISFEDFY QGITAIRNGD PDGQCYGGVA SAQDEEPLAC PDEFDDFVTY EANEVTD SAY MGSESTYSEC ETFTDEDTST LVHPELQPEG DADSAGGS AV PSECLDAMEE PDHGALLLLP GRPHPHGQSV ITVIGGEEHF EDYGEGSEAE LSPETLCNGQ LGCSDPAFLT PSPTKRLSSK KVARYLHQSG ALTMEALEDP SPELMEGPEE DIADKVV FLE RRVLELEKDT AATGEQHSRL RQENLQLVHR ANALEEQLKE QELRACEMVL EETR RQKELL CKMEREKSIE IENLQTRLQQ LDEENSELRS CTPCLKANIE RLEEEKQKLL DEIESLTLRL SEEQENKRRM GDRLSHERHQ FQRDKEATQE LIEDLRKQLE HLQLLKLEAE QRRGRSSSMG LQEYHSRARE SELEQEV RRL KQDNRNLKEQ NEELNGQIIT LSIQGAKSLF STAFSESLAA EISSVSRDEL MEAIQKQEEI NFRLQDYIDR IIVAIMETNP SILEVK
Source	Baculovirus
Target Names	RAB11FIP3
Protein Names	Recommended name: Rab11 family-interacting protein 3 Short name= FIP3- Rab11 Short name= Rab11-FIP3 Alternative name(s): Arfophilin-1 EF hands- containing Rab-interacting protein Short name= Eferin MU-MB-17.148
Expression Region	1-756
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.



Shelf Life

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.

Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.