



Recombinant Human NADH dehydrogenase [ubiquinone] iron-sulfur protein 3, mitochondrial (NDUFS3)

Product Code	CSB-EP015662HU-B
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	O75489
Product Type	Recombinant Protein
Immunogen Species	Homo sapiens (Human)
Purity	≥85% (SDS-PAGE)
Sequence	ESAG ADTRPTVRPR NDVAHKQLSA FGEYVAEILP KYVQQVQVSC FNELEVCIHP DGVIPVLTLF RDHTNAQFKS LVDLTAVDVP TRQNRFEIVY NLLSLRFNSR IRVKTYTDEL TPIESAVSVF KAAWYEREI WDMFGVFFAN HPDLRRILTD YGFEGHPFRK DFPLSGYVEL RYDDEVKRVV AEPVELAQEF RKFDLNSPWE AFPVYRQPPE SLKLEAGDKK PDAK
Source	E.coli
Target Names	NDUFS3
Protein Names	Recommended name: NADH dehydrogenase [ubiquinone] iron-sulfur protein 3, mitochondrial EC= 1.6.5.3 EC= 1.6.99.3 Alternative name(s): Complex I-30kD Short name= CI-30kD NADH-ubiquinone oxidoreductase 30 kDa subunit
Expression Region	37-264
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Full Length of Mature Protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.