



Recombinant Human Zinc finger protein 821 (ZNF821)

Product Code	CSB-MP027066HU
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	O75541
Product Type	Recombinant Protein
Immunogen Species	Homo sapiens (Human)
Purity	>85% (SDS-PAGE)
Sequence	MSRRKQTNPN KVHWDQVFAG LEEQARQAMM KTDFPGDLGS QRQAIQQLRD QDSSSSDSEG DEEETTQDEV SSHTSEEDGG VVKVEKELEN TEQPVGGNEV VEHEVTGNLN SDPLLELCQC PLCQLDCGSR EQLIAHVYQH TAAVVSAKSY MCPVCGRALS SPGSLGRHLL IHSEDQRSNC AVCGARFTSH ATFNSEKLPE VLNMESLPTV HNEGSSAEG KDIAFSPPVY PAGILLVCNN CAAYRKLLEA QTPSVRKWAL RRQNEPLEVR LQRLERERTA KKSRRDNETP EEREVRRMRD REAKRLQRMQ ETDEQRARRL QRDREAMRLK RANETPEKRQ ARLIREREAQ RLKRRLEKMD MMLRAQFGQD PSAMAALAAE MNFFQLPVSG VELDSQLLGK MAFEEQNSSS LH
Source	Mammalian cell
Target Names	ZNF821
Protein Names	Recommended name: Zinc finger protein 821
Expression Region	1-412
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.