



# Recombinant Human Zinc finger protein 821 (ZNF821)

<b>Product Code</b>	CSB-YP027066HU
<b>Storage</b>	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
<b>Uniprot No.</b>	O75541
<b>Product Type</b>	Recombinant Protein
<b>Immunogen Species</b>	Homo sapiens (Human)
<b>Purity</b>	>85% (SDS-PAGE)
<b>Sequence</b>	MSRRKQTNPN KVHWDQVFAG LEEQARQAMM KTDFPGDLGS QRQAIQQLRD QDSSSSDSEG DEEETTQDEV SSHTSEEDGG VVKVEKELEN TEQPVGGNEV VEHEVTGNLN SDPLLELCQC PLCQLDCGSR EQLIAHVYQH TAAVVSAKSY MCPVCGRALS SPGSLGRHLL IHSEDQRSNC AVCGARFTSH ATFNSEKLPE VLNMESLPTV HNEGSSAEG KDIAFSPPVY PAGILLVCNN CAAYRKLLEA QTPSVRKWAL RRQNEPLEVR LQRLERERTA KKSRRDNETP EEREVRRMRD REAKRLQRMQ ETDEQRARRL QRDREAMRLK RANETPEKRQ ARLIREREAQ RLKRRLEKMD MMLRAQFGQD PSAMAALAAE MNFFQLPVSG VELDSQLLGK MAFEEQNSSS LH
<b>Source</b>	Yeast
<b>Target Names</b>	ZNF821
<b>Protein Names</b>	Recommended name: Zinc finger protein 821
<b>Expression Region</b>	1-412
<b>Notes</b>	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
<b>Tag Info</b>	Tag type will be determined during the manufacturing process.
<b>Protein Length</b>	Full length protein
<b>Reconstitution</b>	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
<b>Shelf Life</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.