



# Recombinant Human Forkhead box protein H1 (FOXH1)

<b>Product Code</b>	CSB-BP008817HU
<b>Storage</b>	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
<b>Uniprot No.</b>	O75593
<b>Product Type</b>	Recombinant Protein
<b>Immunogen Species</b>	Homo sapiens (Human)
<b>Purity</b>	>85% (SDS-PAGE)
<b>Sequence</b>	MGPCSGSRLG PPEAESPSQP PKRRKKRYLR HDKPPYTYLA MIALVIQAAP SRRLKLAQII RQVQAVFPFF REDYEGWKDS IRHNLSSNRC FRKVPKDKPAK PQAAGNFWAV DVSLIPAEAL RLQNTALCRR WQNGGARGAF AKDLGPIYVLH GRPYRPPSPP PPPSEGFSIK SLLGGSGEGA PWPGLAPQSS PVPAGTGNSG EEAVPTPPLP SSERPLWPLC PLPGPTRVEG ETVQGGGAI GP STLSPEPRAW PLHLLQGTA V PGGRSSGGHR ASLWGQLPTS YLPIYTPNVV MPLAPPPTSC PQCSTSPAY WGVAPETRGP PGLLCDLDAL FQGVPPNKSI YDVVWVSHPRD LAAPGPGWLL SWCSL
<b>Source</b>	Baculovirus
<b>Target Names</b>	FOXH1
<b>Protein Names</b>	Recommended name: Forkhead box protein H1 Alternative name(s): Forkhead activin signal transducer 1 Short name= Fast-1 Short name= hFAST-1 Forkhead activin signal transducer 2 Short name= Fast-2
<b>Expression Region</b>	1-365
<b>Notes</b>	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
<b>Tag Info</b>	Tag type will be determined during the manufacturing process.
<b>Protein Length</b>	Full length protein
<b>Target Details</b>	FOXH1 encodes a human homolog of Xenopus forkhead activin signal transducer-1. FOXH1 protein binds SMAD2 and activates an activin response element via binding the DNA motif TGT(G/T)(T/G)ATT.
<b>Reconstitution</b>	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
<b>Shelf Life</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.