



Recombinant Human Trafficking protein particle complex subunit 6A (TRAPPC6A)

Product Code	CSB-MP024232HU
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	O75865
Product Type	Recombinant Protein
Immunogen Species	Homo sapiens (Human)
Purity	≥85% (SDS-PAGE)
Sequence	MADTVLFEFL HTEMVAELWA HDPDPGPGGQ KMSLSVLEGM GFRVGQALGE RLPRETLAFR EELDVLKFLC KDLWVAVFQK QMDSLRTNHQ GTYVLQDNSF PLLLPASGL QYLEEAPKFL AFTCGLLRGA LYTLGIESVV TASVAALPVC KFQVVIPKS
Source	Mammalian cell
Target Names	TRAPPC6A
Protein Names	Recommended name: Trafficking protein particle complex subunit 6A Short name= TRAPP complex subunit 6A
Expression Region	1-159
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.