



Recombinant Bovine Neutrophil cytosol factor 1 (NCF1)

Product Code	CSB-YP015525BO
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	O77774
Product Type	Recombinant Protein
Immunogen Species	Bos taurus (Bovine)
Purity	>85% (SDS-PAGE)
Sequence	MGDHFIRHIA LLGFEKRFVP SQHYVYMFLV KWQDLSEKVV YRRFTEIYEF HKILKEMFPI EAGDINPENR IIPHLPAWRW YDGQRVAESR QGTLTEYCST LMSLPVKISR CPHLLNFFKV RPDDLKLPD SQVKKPETYL MPKDGKNNAA DITGPILQT YRAIADYEKG SSSQMALATG DVVDVVEKNE SGWWFCQMKT KRGWVPASYL EPLDSPDEAE DPEPNYAGEP YITIKAYTAV LEDEISLEEG EAIEVIHKLK DGWWVIRKED VTGYFPSMYL QKAGQDVAQA KSQIKSRGAP PRRSSIRNAH SIHQRSRKRL SQDTYRRNSV RFMQQRRHQ R PEPQRSRSAL REQQQPKTER PKPQPAVPPR PSADLILHRC SESTRKRLAS AV
Source	Yeast
Target Names	NCF1
Protein Names	Recommended name: Neutrophil cytosol factor 1 Short name= NCF-1 Alternative name(s): 47 kDa neutrophil oxidase factor NCF-47K Neutrophil NADPH oxidase factor 1 p47-phox
Expression Region	1-392
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.