



Recombinant Bovine Angiopoietin-2 (ANGPT2)

Product Code	CSB-MP001707BO
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	O77802
Product Type	Recombinant Protein
Immunogen Species	Bos taurus (Bovine)
Purity	≥85% (SDS-PAGE)
Sequence	HS GSRKGMIDIAA GKKQYQVQHG ACSYTFLLPE TDHCRSPSSA YVPNAVQRDA PLDYDDSVQR LQVLENIMEN NTQWLMKLEN YIQDNMKKEM VEIQQNAVQN QTAVMIEIGT NLLNQTAEQT RKLTDVEAQV LNQTTRLELQ LLEHSLSTNK LEKQILDQTS EISKLDKNS FLEKKVLDME DKHIVQLRSI KEEKDQLQVL VSKQNSIIEE LEKQLVTATV NNSVLQKQQH DLMETVNNLL TLMSTSNPSY SLLAKDEQII FRDCGEAFKS GLTTSGVYTL TFPNSTEEIK AYCDMETGGG GWTVIQRRED GSVDFQRTWK EYKVGFGNPS GEHWLGNEFV SQVTGQKRYV LKIHRLDWEH NEAYSLEYDHF YLSNEELNYR IHLKGLTGTA GKISSISQPG NDFSTKDADN DKCICKCSQM LTGGWWFDAC GPSNLNGMY YPQRQNTNKF N GIKWYYWKGS GYSLKATTMM IRPADF
Source	Mammalian cell
Target Names	ANGPT2
Protein Names	Recommended name: Angiopoietin-2 Short name= ANG-2
Expression Region	19-496
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Full Length of Mature Protein
Target Details	This protein is an antagonist of angiopoietin 1 (ANGPT1) and endothelial TEK tyrosine kinase (TIE-2, TEK). The encoded protein disrupts the vascular remodeling ability of ANGPT1 and may induce endothelial cell apoptosis. Three transcript variants encoding three different isoforms have been found for this gene.
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.