



Recombinant Mouse Peptidoglycan recognition protein 1 (Pglyrp1)

Product Code	CSB-EP017862MO-B
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	O88593
Product Type	Recombinant Protein
Immunogen Species	Mus musculus (Mouse)
Purity	>85% (SDS-PAGE)
Sequence	FI VPRSEWRALP SECSSRLGHP VRYVVISHTA GSFCNSPDSC EQQARNVQHY HKNELGWCDV AYNFLIGEDG HVYEGRGWNI KGDHTGPIWN PMSIGITFMG NFMDRVPKR ALRAALNLE CGVSRGFLRS NYEVKGHRDV QSTLSPGDQL YQVIQSWEHY RE
Source	E.coli
Target Names	Pglyrp1
Protein Names	Recommended name: Peptidoglycan recognition protein 1 Alternative name(s): Cytokine tag7 Peptidoglycan recognition protein short Short name= PGRP-S
Expression Region	19-182
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Full Length of Mature Protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.