



Recombinant Mouse E3 SUMO-protein ligase PIAS1 (Pias1), partial

Product Code	CSB-BP017957MO1
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	O88907
Product Type	Recombinant Protein
Immunogen Species	Mus musculus (Mouse)
Purity	>85% (SDS-PAGE)
Sequence	MADSAELKQMVMVSLRVSELQVLLGYAGR NKHGRKHELLTKALHLLKAGCSPA VQMKIKELYRRRFPQKIMTPADLSIPNVHSSPMPPTLSPSTIPQLTYDGHPASS PLLPSVLLGPKHELELPHLTSALHPVHPDIK LQKLPFYDLLDELIKPTSLASDNS QRFRETCFAFALTPQQVQQISSSMDISG TKCDFTVQVQLRFCLSETSCPQEDH FPPNLCVKVNTKPCSLPGYLPPTKNGVEPKRPSRPINITSLVRLSTTVPNTIVVS WTAEIGRTYSMAVYLVKQLSSTVLLQRLRAKGIRNPDHSRALIKEKLTADPDSEI ATTSLRVSLLCPLGKMRLTIPCRALTC SHLQCFDATLYIQMNEKKPTWVCPVC DKKAPYEHLIIDGLFMEILKYCT
Source	Baculovirus
Target Names	Pias1
Protein Names	Recommended name: E3 SUMO-protein ligase PIAS1 EC= 6.3.2.-Alternative name(s): DEAD/H box-binding protein 1 Protein inhibitor of activated STAT protein 1
Expression Region	1-400
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Partial
Target Details	This gene encodes a member of the mammalian PIAS [protein inhibitor of activated STAT-1 (signal transducer and activator of transcription-1)] family. This member contains a putative zinc-binding motif and a highly acidic region. It inhibits STAT1-mediated gene activation and the DNA binding activity, binds to Gu protein/RNA helicase II/DEAD box polypeptide 21, and interacts with androgen receptor (AR). It functions in testis as a nuclear receptor transcriptional coregulator and may have a role in AR initiation and maintenance of spermatogenesis.
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.



Shelf Life

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.

Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.