



Recombinant Mouse C-C motif chemokine 20 (Ccl20)

Product Code	CSB-BP004784MO
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	O89093
Product Type	Recombinant Protein
Immunogen Species	Mus musculus (Mouse)
Purity	>85% (SDS-PAGE)
Sequence	ASN YDCCLSYIQT PLPSRAIVGF TRQMADEACD INAIIFHTKK RKSVCADPKQ NWWKRAVNLL SLRVKKM
Source	Baculovirus
Target Names	Ccl20
Protein Names	Recommended name: C-C motif chemokine 20 Alternative name(s): Beta-chemokine exodus-1 CC chemokine LARC CC chemokine ST38 Liver and activation-regulated chemokine Macrophage inflammatory protein 3 alpha Short name= MIP-3- α p
Expression Region	28-97
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Full Length of Mature Protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.