



Recombinant Human HAUS augmin-like complex subunit 5 (HAUS5)

Product Code	CSB-MP012212HU
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	O94927
Product Type	Recombinant Protein
Immunogen Species	Homo sapiens (Human)
Purity	≥85% (SDS-PAGE)
Sequence	MELAQEAREL GCWAVEEMGV PVAARAPEST LRRCLGQGA DIWAYILQHV HSQRTVKKIR GNLLWYGHQD SPQVRRKLEL EAAVTRLRAE IQELDQSLEL MERDTEAQDT AMEQARQHTQ DTQRRALLR AQAGAMRRQQ HTLRDPMQRL QNQLRRLQDM ERKAKVDVTF GSLSAALGL EPVVLRDVRT ACTLRAQFLQ NLLLPQAKRG SLPTPHDDHF GTSYQQWLSS VETLLTNHPP GHVLALEHL AAERAEIRS LCSGDGLGDT EISRPQAPDQ SDSSQTLPSM VHLIQEGWRT VGVLVSRST LLKERQVLTQ RLQGLVEEVE RRVLGSSERQ VLILGLRRCC LWTELKALHD QSQELQDAAG HRQLLLRELQ AKQQRILHWR QLVEETQEQQV RLLIKGNSAS KTRLCRSPGE VLALVQRKVV PTFEAVAPQS RELLRCLEEE VRHLPHILLG TLLRHRPGEL KPLPTVLP SI HQLHPASPRG SSFIALSHKL GLPPGKASEL LLPAAASLRQ DLLLLQDQRS LWCWDLHMK TSLPPGLPTQ ELLQIQASQE KQKQENLGQA LKRLEKLLKQ ALERIPELQG IVGDWWEQPG QAALSEELCQ GLSLPQWRLR WVQAQGALQK LCS
Source	Mammalian cell
Target Names	HAUS5
Protein Names	Recommended name: HAUS augmin-like complex subunit 5
Expression Region	1-633
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.