



# Recombinant Human HAUS augmin-like complex subunit 5 (HAUS5)

|                          |                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                               |
|--------------------------|---------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| <b>Product Code</b>      | CSB-YP012212HU                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                |
| <b>Storage</b>           | Store at -20°C, for extended storage, conserve at -20°C or -80°C.                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                             |
| <b>Uniprot No.</b>       | O94927                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                        |
| <b>Product Type</b>      | Recombinant Protein                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                           |
| <b>Immunogen Species</b> | Homo sapiens (Human)                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                          |
| <b>Purity</b>            | ≥85% (SDS-PAGE)                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                               |
| <b>Sequence</b>          | <p>MELAQEAREL GCWAVEEMGV PVAARAPEST LRRCLGQGA DIWAYILQHV<br/>           HSQRTVKKIR GNLLWYGHQD SPQVRRKLEL EAAVTRLRAE IQELDQSLEL<br/>           MERDTEAQDT AMEQARQHTQ DTQRRALLR AQAGAMRRQQ<br/>           HTLRDPMQRL QNQLRRLQDM ERKAKVDVTF GSLSAALGL EPVVLRDVRT<br/>           ACTLRAQFLQ NLLLPQAKRG SLPTPHDDHF GTSYQQWLSS VETLLTNHPP<br/>           GHVLALEHL AAERAEIRS LCSGDGLGDT EISRPQAPDQ SDSSQTLPSM<br/>           VHLIQEGWRT VGVLVSRST LLKERQVLTQ RLQGLVEEVE RRVLGSSERQ<br/>           VLILGLRRCC LWTELKALHD QSQELQDAAG HRQLLLRELQ AKQQRILHWR<br/>           QLVEETQEQV RLLIKGNSAS KTRLCRSPGE VLALVQRKVV PTFEAVAPQS<br/>           RELLRCLEEE VRHLPHILLG TLLRHRPGEL KPLPTVLP SI HQLHPASPRG<br/>           SSFIALSHKL GLPPGKASEL LLPAAASLRQ DLLLLQDQRS LWCWDLHMK<br/>           TSLPPGLPTQ ELLQIQASQE KQKQENLGQA LKRLEKLLKQ ALERIPELQG<br/>           IVGDWWEQPG QAALSEELCQ GLSLPQWRLR WVQAQGALQK LCS</p> |
| <b>Source</b>            | Yeast                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                         |
| <b>Target Names</b>      | HAUS5                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                         |
| <b>Protein Names</b>     | Recommended name: HAUS augmin-like complex subunit 5                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                          |
| <b>Expression Region</b> | 1-633                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                         |
| <b>Notes</b>             | Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                           |
| <b>Tag Info</b>          | Tag type will be determined during the manufacturing process.                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                 |
| <b>Protein Length</b>    | Full length protein                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                           |
| <b>Reconstitution</b>    | We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                   |
| <b>Shelf Life</b>        | <p>The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.</p> <p>Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.</p>                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                 |