



Recombinant Human Aflatoxin B1 aldehyde reductase member 3 (AKR7A3)

Product Code	CSB-MP001551HU
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	O95154
Product Type	Recombinant Protein
Immunogen Species	Homo sapiens (Human)
Purity	≥85% (SDS-PAGE)
Sequence	MSRQLSRARP ATVLGAMEMG RRMDAPTSAA VTRAFLEIGH TEIDTAFVYS EGQSETILGG LGLRLGGSDC RVKIDTKAIP LFGNSLKPDS LRFQLETSLK RLQCPRVDLF YLHMPDHSTP VEETLRACHQ LHQEGKFVEL GLSNYAWEV AEICTLCKSN GWILPTVYQG MYNAITRQVE TELFPCLRHF GLRFYAFNPL AGLLTGKYK YEDKNGKQPV GRFFGNTWAE MYRNRWKEH HFEGIALVEK ALQAAYGASA PSMTSATLRW MYHHSQ LQGA HGDAVILGMS SLEQLEQNLA AAEEGPLIPA VVDAFNQAWH LVTHECPNYF R
Source	Mammalian cell
Target Names	AKR7A3
Protein Names	Recommended name: Aflatoxin B1 aldehyde reductase member 3 EC= 1.-.-.- Alternative name(s): AFB1 aldehyde reductase 2 Short name= AFB1-AR 2
Expression Region	1-331
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.