



Recombinant Bovine Inositol-tetrakisphosphate 1-kinase (ITPK1)

Product Code	CSB-MP011908BO
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	P0C0T1
Product Type	Recombinant Protein
Immunogen Species	Bos taurus (Bovine)
Purity	≥85% (SDS-PAGE)
Sequence	<p>MQTFLKGKRV GYWLSEKKIK KLNFAQFAEL CRKRGIEVVQ LNLSRPIEEQ GPLDVIIHKL TDVILEADQN DSQALELVHR FQEYIDAHPE TIVLDPLPAI RTLLDRSKSY ELIRKIEAYM KDDRICSPPF MELTSLCGDD TMRLLLEENGL AFPFICKTRV AHGTNSHEMA IVFNQEGLSA IQPPCVVQNF INHNAVLYKV FVVGESYTVV QRPSLKNFSA GTSDRESIFF NSHNVSKPES SSVLTALDKI EGVFERPSDE VIRELSRALR QALGVSLFGI DIIINNQTGQ HAVIDINAFP GYEGVSEFFT DLLNHASVL QGQSSGVAGA GDVAPLKHSR LLAEQAGGLA AERTCSASPG CCSSMMGQEP PWTPEADMGG VGAGSTAKLP HQRLGCTAGV SPSFQQHCVA SLATKASSQ</p>
Source	Mammalian cell
Target Names	ITPK1
Protein Names	<p>Recommended name: Inositol-tetrakisphosphate 1-kinase EC= 2.7.1.134 Alternative name(s): Inositol 1,3,4-trisphosphate 5/6-kinase Short name= Inositol-triphosphate 5/6-kinase Short name= Ins(1,3,4)P(3) 5/6-kinase EC= 2.7</p>
Expression Region	1-419
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	<p>The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.</p> <p>Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.</p>