



Recombinant Human Uridine diphosphate glucose pyrophosphatase (NUDT14)

Product Code	CSB-BP016159HU
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	O95848
Product Type	Recombinant Protein
Immunogen Species	Homo sapiens (Human)
Purity	>85% (SDS-PAGE)
Sequence	MERIEGASVG RCAASPYLRP LTLHYRQNGA QKSWDFMKTH DSVTVLLFNS SRRSLVLVKQ FRPAVYAGEV ERRFPGLAA VDQDGPRELQ PALPGSAGVT VELCAGLVDQ PGLSLEEVAC KEAWEECGYH LAPSDLRRA TYWSGVGLTG SRQTMFYTEV TDAQRSRGG GLVEEGELIE VVHLPLEGAQ AFADDPDIPK TLGVIFGVSF FLSQVAPNLD LQ
Source	Baculovirus
Target Names	NUDT14
Protein Names	Recommended name: Uridine diphosphate glucose pyrophosphatase Short name= UDPG pyrophosphatase Short name= UGPPase EC= 3.6.1.45 Alternative name(s): Nucleoside diphosphate-linked moiety X motif 14 Short name= Nudix moti
Expression Region	1-222
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.