



Recombinant Human Lymphocyte antigen 6 complex locus protein G6c (LY6G6C)

Product Code	CSB-BP013245HU
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	O95867
Product Type	Recombinant Protein
Immunogen Species	Homo sapiens (Human)
Purity	>85% (SDS-PAGE)
Sequence	DI RCHSCYKVPV LGCVDRQSCR LEPGQQCLTT HAYLGKMWVF SNLRCGTPEE PCQEAFNQTN RKLGLTYNTT CCNKDNCNS
Source	Baculovirus
Target Names	LY6G6C
Protein Names	Recommended name: Lymphocyte antigen 6 complex locus protein G6c
Expression Region	19-99
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Full Length of Mature Protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.