



Recombinant Human Pre-mRNA-splicing factor SYF2 (SYF2)

| | |
|--------------------------|---|
| Product Code | CSB-EP022996HU |
| Storage | Store at -20°C, for extended storage, conserve at -20°C or -80°C. |
| Uniprot No. | O95926 |
| Product Type | Recombinant Protein |
| Immunogen Species | Homo sapiens (Human) |
| Purity | >85% (SDS-PAGE) |
| Sequence | AAIAASEVL VDSAEEGSLA AAAELAAQKR EQRLRKFREL HLMRNEARKL NHQEVVEEDK RLKLPANWEA K KARLEWELK EEEKKKECAA RGEDYEKVKL LEISAEDAER WERKKRKNP DLGFSDYAAA QLRQYHRLTK QIKPDMETYE RLREKHGEEF FPTSNSLLHG THVPSTEEID RMVIDLEKQI EKRDKYSRRR PYNDADIDY INERNAKFNK KAERFYGKYT AEIKQNLERG TAV |
| Source | E.coli |
| Target Names | SYF2 |
| Protein Names | Recommended name: Pre-mRNA-splicing factor SYF2 Alternative name(s): CCNDBP1-interactor p29 |
| Expression Region | 2-243 |
| Notes | Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week. |
| Tag Info | Tag type will be determined during the manufacturing process. |
| Protein Length | Full Length of Mature Protein |
| Target Details | This gene encodes a nuclear protein that interacts with cyclin D-type binding-protein 1, which is thought to be a cell cycle regulator at the G1/S transition. Alternate transcriptional splice variants, encoding different isoforms, have been characterized. |
| Reconstitution | We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference. |
| Shelf Life | The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C. |