



# Recombinant Human Pre-mRNA-splicing factor SYF2 (SYF2)

|                          |   |
|--------------------------|---|
| <b>Product Code</b>      | CSB-YP022996HU  |
| <b>Storage</b>           | Store at -20°C, for extended storage, conserve at -20°C or -80°C.   |
| <b>Uniprot No.</b>       | O95926  |
| <b>Product Type</b>      | Recombinant Protein   |
| <b>Immunogen Species</b> | Homo sapiens (Human)  |
| <b>Purity</b>            | ≥85% (SDS-PAGE)   |
| <b>Sequence</b>          | AAIAASEVL VDSAEEGSLA AAAELAAQKR EQRLRKFREL HLMRNEARKL<br>NHQEVVEEDK RLKLPANWEA K KARLEWELK EEEKKKECAA RGEDYEKVKL<br>LEISAEDAER WERKKKRKNP DLGFSDYAAA QLRQYHRLTK QIKPDMETYE<br>RLREKHGEEF FPTSNSLLHG THVPSTEEID RMVIDLEKQI EKRDKYSRRR<br>PYNDADIDY INERNAKFNK KAERFYGKYT AEIKQNLERG TAV  |
| <b>Source</b>            | Yeast   |
| <b>Target Names</b>      | SYF2  |
| <b>Protein Names</b>     | Recommended name: Pre-mRNA-splicing factor SYF2 Alternative name(s):<br>CCNDBP1-interactor p29  |
| <b>Expression Region</b> | 2-243   |
| <b>Notes</b>             | Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.   |
| <b>Tag Info</b>          | Tag type will be determined during the manufacturing process.   |
| <b>Protein Length</b>    | Full Length of Mature Protein   |
| <b>Target Details</b>    | This gene encodes a nuclear protein that interacts with cyclin D-type binding-protein 1, which is thought to be a cell cycle regulator at the G1/S transition. Alternate transcriptional splice variants, encoding different isoforms, have been characterized.   |
| <b>Reconstitution</b>    | We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference. |
| <b>Shelf Life</b>        | The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.   |