



# Recombinant Human BUB3-interacting and GLEBS motif-containing protein ZNF207 (ZNF207)

<b>Product Code</b>	CSB-MP026578HU
<b>Storage</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
<b>Uniprot No.</b>	O43670
<b>Storage Buffer</b>	Lyophilized from Tris/PBS-based buffer, 6% Trehalose, pH 8.0
<b>Product Type</b>	Recombinant Proteins
<b>Immunogen Species</b>	Homo sapiens (Human)
<b>Purity</b>	>85% (SDS-PAGE)
<b>Sequence</b>	MGRKKKKQLK PWCWYCNRDF DDEKILIQHQ KAKHFKCHIC HKKLYTGPGL AIHCMQVHKE TIDAVPNAIP GRTDIELEIY GMEGIPEKDM DERRRLEQK TQESQKKKQQ DDSDEYDDDD SAASTSFQPQ PVQPQQGYIP PMAQPGLPPV PGAPGMPPGI PPLMPGVVPL MPGMPPVMPG MPPGMMPMGG MMPPGPGIPP LMPGMPPGMP PPVPRPGIPP MTQAQAVSAP GILNRPPAPT ATVPAQPPV TKPLFPSAGQ MGTPVTSSST ASSNSELSA SSKALFPSTA QAQAAVQGPV GTDFKPLNST PATTTEPPKP TFPAYTQSTA STTSTNSTA AKPAASITSK PATLTTTSAT SKLIHPDEDI SLEERRAQLP KYQRNLPRPG QAPIGNPPVG PIGGMMPPQP GIPQQQGMRP PMPPHGQYGG HHQGMPPGYP GAMPPYQGGP PMVPPYQGGP PRPPMGMRPP VMSQGGRY
<b>Source</b>	Mammalian cell
<b>Target Names</b>	ZNF207
<b>Protein Names</b>	Recommended name: Zinc finger protein 207
<b>Expression Region</b>	1-478
<b>Notes</b>	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
<b>Tag Info</b>	Tag type will be determined during the manufacturing process.
<b>Protein Length</b>	Full length protein
<b>Reconstitution</b>	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
<b>Shelf Life</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life



of lyophilized form is 12 months at -20°C/-80°C.