



Recombinant Human BUB3-interacting and GLEBS motif-containing protein ZNF207 (ZNF207)

Product Code	CSB-MP026578HU
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	O43670
Storage Buffer	Lyophilized from Tris/PBS-based buffer, 6% Trehalose, pH 8.0
Product Type	Recombinant Proteins
Immunogen Species	Homo sapiens (Human)
Purity	≥85% (SDS-PAGE)
Sequence	MGRKKKKQLK PWCWYCNRDF DDEKILIQHQ KAKHFKCHIC HKKLYTGPGL AIHCMQVHKE TIDAVPNAIP GRDIELEIY GMEGIPEKDM DERRRLEEQK TQESQKKKQQ DDSDEYDDDD SAASTSFQPQ PVQPQQGYIP PMAQPGLPPV PGAPGMPPGI PPLMPGVPL MPGMPVMPG MPPGMMPMGG MMPPGPGIPP LMPGMPPGMP PPVPRGIPP MTQAQAVSAP GILNRPPAPT ATVPAQPPV TKPLFPSAGQ MGTPVTSSST ASSNSELSA SSKALFPSTA QAQAAVQGPV GTDFKPLNST PATTTEPPKP TFPAYTQSTA STTSTNSTA AKPAASITSK PATLTTTSAT SKLIHPDEDI SLEERRAQLP KYQRNLPRPG QAPIGNPPVG PIGGMMPPQP GIPQQGMRP PMPPHGQYGG HHQGMPGYLP GAMPPYGQGP PMVPPYQGGP PRPPMGMRPP VMSQGGRY
Source	Mammalian cell
Target Names	ZNF207
Protein Names	Recommended name: Zinc finger protein 207
Expression Region	1-478
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life



of lyophilized form is 12 months at -20°C/-80°C.