



Recombinant Human U3 small nucleolar RNA-interacting protein 2 (RRP9)

Product Code	CSB-BP020532HU
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	O43818
Product Type	Recombinant Protein
Immunogen Species	Homo sapiens (Human)
Purity	≥85% (SDS-PAGE)
Sequence	MSATAAARKR GKPASGAGAG AGAGKRRRKA DSAGDRGKSK GGGKMNEEIS SDSESESLAP RKPEEEEEEE LEETAQEKKL RLAKLYLEQL RQQEEEEKAEA RAFEEEDQVAG RLKEDVLEQR GRLQKLVAKE IQAPASADIR VLRGHQLSIT CLVVTTPDDSA IFSAAKDCSI IKWSVESGRK LHVIPRAKKG AEGKPPGHSS HVLCAISSD GKYLASGDRS KLILIWEAQS CQHLYTFTGH RDAVSGLAFR RGTHQLYSTS HDRSVKVWNV AENSYVETLF GHQDAVAALD ALSRECCVTA GGRDGTVRVW KIPEESQLVF YGHQGSIDCI HLINEEHMVS GADDGSVALW GLSKKRPLAL QREAHGLRGE PGLEQPFWIS SVAALLNTDL VATGSHSSCV RLWQCGEGFR QLDLLCDIPL VGFINSLKFS SSGDFLVAGV GQEHRLGRWW RIKEARN SVC IIPLRVPVP PAAGS
Source	Baculovirus
Target Names	RRP9
Protein Names	Recommended name: U3 small nucleolar RNA-interacting protein 2 Alternative name(s): RRP9 homolog U3 small nucleolar ribonucleoprotein-associated 55 kDa protein Short name= U3 snoRNP-associated 55 kDa protein Short name= U3-55K
Expression Region	1-475
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.