



# Recombinant Human Beta-1,4-galactosyltransferase 3 (B4GALT3), partial, Yeast

<b>Product Code</b>	CSB-YP002515HU
<b>Storage</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
<b>Uniprot No.</b>	O60512
<b>Product Type</b>	Recombinant Protein
<b>Immunogen Species</b>	Homo sapiens (Human)
<b>Purity</b>	≥85% (SDS-PAGE)
<b>Source</b>	Yeast
<b>Target Names</b>	B4GALT3
<b>Protein Names</b>	Recommended name: Beta-1,4-galactosyltransferase 3 Short name= Beta-1,4-GalTase 3 Short name= Beta4Gal-T3 Short name= b4Gal-T3 EC= 2.4.1.- Alternative name(s): UDP-Gal:beta-GlcNAc beta-1,4-galactosyltransferase 3 UDP-galact
<b>Notes</b>	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
<b>Tag Info</b>	N-terminal His-tagged Tag-Free
<b>Protein Length</b>	Full length protein
<b>Shelf Life</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.