



Recombinant Human Heat shock factor 2-binding protein (HSF2BP)

Product Code	CSB-MP010793HU
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	O75031
Product Type	Recombinant Protein
Immunogen Species	Homo sapiens (Human)
Purity	≥85% (SDS-PAGE)
Sequence	MGEAGAAEEA CRHMGTKKEF VKVRKKDLER LTTEVMQIRD FLPRILNGEV LESFQKLKIV EKNLERKEQE LEQLKMDCEH FKARLETVQA DNIREKKEKL ALRQQLNEAK QQLLQQAIEYC TEMGAACTL LWGVSSSEEV VKAILGGDKA LKFFSITGQT MESFVKSLDG DVQELDSDES QFVFALAGIV TNVAAIACGR EFLVNSSRVL LDTILQLLGD LKPGQCTKLK VLMLMSLYNV SINLKGLKYI SESPGFIPLL WWLLSDPDAE VCLHVLRLVQ SVVLEPEVFS KSASEFRSSL PLQRILAMSK SRNPRLQTAA QELLEDLRTL EHNV
Source	Mammalian cell
Target Names	HSF2BP
Protein Names	Recommended name: Heat shock factor 2-binding protein
Expression Region	1-334
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Full length protein
Target Details	HSF2 binding protein (HSF2BP) associates with HSF2. The interaction occurs between the trimerization domain of HSF2 and the amino terminal hydrophilic region of HSF2BP that comprises two leucine zipper motifs. HSF2BP may therefore be involved in modulating HSF2 activation.
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.