



Recombinant Human KH domain-containing, RNA-binding, signal transduction-associated protein 3 (KHDRBS3)

Product Code	CSB-MP012156HU
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	O75525
Product Type	Recombinant Protein
Immunogen Species	Homo sapiens (Human)
Purity	>85% (SDS-PAGE)
Sequence	MEEKYLPPELM AEKDSLDPFSF THALRLVNQE IEKFQKGEGK DEEKYIDVVI NKNMKLGQKV LIPVKQFPKF NFVGKLLGPR GNSLKRLQEE TLTKMSILGK GSMRDKAKEE ELRKSGEAKY FHLNDDLHVL IEVFAPPAEA YARMGHAAEE IKKFLIPDYN DEIRQAQLQE LTYLNGGSEN ADVPVVRGKP TLTRGVVPAP AITRGRGGVT ARPVGVVVPR GTPTPRGVLS TRGPVSRGRG LLTPRARGVP PTGYRPPPPP PTQETYGEYD YDDGYGTAYD EQSYDSYDNS YSTPAQSGAD YYDYGHGLSE ETYDSYGQEE WTNSRHKAPS ARTAKGVYRD QPYGRY
Source	Mammalian cell
Target Names	KHDRBS3
Protein Names	Recommended name: KH domain-containing, RNA-binding, signal transduction-associated protein 3 Alternative name(s): RNA-binding protein T-Star Sam68-like mammalian protein 2 Short name= SLM-2 Sam68-like phosphotyrosine protein
Expression Region	1-346
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.