



Recombinant Human E3 SUMO-protein ligase PIAS1 (PIAS1)

Product Code	CSB-EP017957HU-B
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	O75925
Product Type	Recombinant Protein
Immunogen Species	Homo sapiens (Human)
Purity	>85% (SDS-PAGE)
Sequence	<p>ADSAELKQM VMSLRVSELQ VLLGYAGRNK HGRKHELLTK ALHLLKAGCS PAVQMKIKEL YRRRFPQKIM TPADLSIPNV HSSPMPATLS PSTIPQLTYD GHPASSPLLP VSLLGPKHEL ELPHLTSALH PVHPDIKLQK LPFYDLLDEL IKPTSLASDN SQRFRETCA FALTPQQVQQ ISSSMDISGT KCDFTVQVQL RFCLSETSCP QEDHFPPNLC VKVNTKPCSL PGYLPPTKNG VEPKRPSRPI NITSLVRLST TVPNTIVVSW TAEIGRNYSM AVYLVKQLSS TVLLQRLRAK GIRNPDHSRA LIKEKLTADP DSEIATTSR VSLLCPLGKM RTIPCRALT CSHLQCFDAT LYIQMNEKKP TWVCPVCDKK APYEHLIDG LFMEILKYCT DCDEIQFKED GTWAPMRSKK EVQEVSASYN GVDGCLSSTL EHQVASHHQS SNKNKKVEVI DLTIDSSSDE EEEEPSAKRT CPSLSPTSPL NNGILSLPH QASPVSRTPS LPAVDTSYIN TSLIQDYRHP FHMTMPYDL QGLDFFPFLS GDNQHYNTSL LAAAAAAVSD DQDLLHSSRF FPYTSSQMFL DQLSAGGSTS LPTTNGSSSG SNSSLVSSNS LRESHSHVT NRSSTDTASI FGIIPDIISL D</p>
Source	E.coli
Target Names	PIAS1
Protein Names	Recommended name: E3 SUMO-protein ligase PIAS1 EC= 6.3.2.- Alternative name(s): DEAD/H box-binding protein 1 Gu-binding protein Short name= GBP Protein inhibitor of activated STAT protein 1 RNA helicase II-binding protein
Expression Region	2-651
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Full Length of Mature Protein
Target Details	<p>This gene encodes a member of the mammalian PIAS [protein inhibitor of activated STAT-1 (signal transducer and activator of transcription-1)] family. This member contains a putative zinc-binding motif and a highly acidic region. It inhibits STAT1-mediated gene activation and the DNA binding activity, binds to Gu protein/RNA helicase II/DEAD box polypeptide 21, and interacts with androgen receptor (AR). It functions in testis as a nuclear receptor transcriptional coregulator and may have a role in AR initiation and maintenance of spermatogenesis.</p>



Reconstitution

We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.

Shelf Life

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.