



Recombinant Human Serine/threonine-protein kinase 17B (STK17B)

Product Code	CSB-MP022838HU
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	O94768
Product Type	Recombinant Protein
Immunogen Species	Homo sapiens (Human)
Purity	>85% (SDS-PAGE)
Sequence	MSRRRFDCRS ISGLLTTPQ IPIKMENFNN FYILTSKELG RGKFAVVRQC ISKSTGQEYA AKFLKKRRRG QDCRAEILHE IAVLELAKSC PRVINLHEVY ENTSEIILIL EYAAGGEIFS LCLPELAEMV SENDVIRLIK QILEGVYYLH QNNIVHLDLK PQNILLSSIY PLGDIKIVDF GMSRKIGHAC ELREIMGTP YLAPEILNYD PITTATDMWN IGIAYMLLT HTSPFVGEDN QETYLNISQV NVDYSEETFS SVSQLATDFI QSLLVKNPEK RPTAEICLSH SWLQQWDFEN LFHPEETSSS SQTQDHSVRS SEDKTSKSSC NGTCGDREDK ENIPEDSSMV SKRFRFDDSL PNPHELVS DL LC
Source	Mammalian cell
Target Names	STK17B
Protein Names	Recommended name: Serine/threonine-protein kinase 17B EC= 2.7.11.1 Alternative name(s): DAP kinase-related apoptosis-inducing protein kinase 2
Expression Region	1-372
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.