



# Recombinant Human Dendrin (DDN)

<b>Product Code</b>	CSB-MP006592HU
<b>Storage</b>	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
<b>Uniprot No.</b>	O94850
<b>Product Type</b>	Recombinant Protein
<b>Immunogen Species</b>	Homo sapiens (Human)
<b>Purity</b>	≥85% (SDS-PAGE)
<b>Sequence</b>	MLDGPLFSEG PDSPRELQDE ESGSCLWVQK SKLLVIEVKT ISCHYSRRAP SRQPMDFQAS HWARGFQNRG CGPRPGSPQP PRRRPWASRV LQEATNWRAG PLAEVAREQ EKRKAASQER EAKETERKRR KAGGARRSPP GRPRPEPRNA PRVAQLAGLP APLRPERLAP VGRAPRPSAQ PQSDPGSAWA GPWGGRRPGP PSYEAHLLLR GSAGTAPRRR WDRPPPYVAP PSYEGPHRTL GTKRGPNSQ VPTSSAPAAT PARTDGGRTK KRLDPRIYRD VLGAWGLRQG QGLLGGSPGC GAARARPEPG KGVVEKSLGL AAADLNSGSD SHPQAKATGS AGTEIAPAGS ATAAPCAPHP APRSRHHLKG SREGKEGEQI WFPKCWIPSP KKQPPRHSQT LPRPWAPGGT GWRESLGLGE GAGPETLEGW KATRRRAHTLP RSSQGLSRGE GVFVIDATCV VIRSQYVPTP RTQQVQLLPS GVTRVVGDSP SQSKPGKEEG EGATVFPSPC QKRLSSSRLH HQPGGGRGGE AEGGRPGDST LEERTFRILG LPAPEVNLRD APTQPGSPEH QALGPAASGA QGRAEGSEVA VVQRRAGRGW ARTPGPYAGA LREAVSRIRR HTAPDSDTDE AEELSVHSGS SDGSDTEAPG ASWRNERTLP EVGNSSPEED GKTAELSDSV GEILDVISQT EEVLFGVRDI RGTQQGNRKR Q
<b>Source</b>	Mammalian cell
<b>Target Names</b>	DDN
<b>Protein Names</b>	Recommended name: Dendrin
<b>Expression Region</b>	1-711
<b>Notes</b>	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
<b>Tag Info</b>	Tag type will be determined during the manufacturing process.
<b>Protein Length</b>	Full length protein
<b>Reconstitution</b>	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
<b>Shelf Life</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.