



# Recombinant Human Zinc finger protein 432 (ZNF432)

<b>Product Code</b>	CSB-EP026746HU
<b>Storage</b>	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
<b>Uniprot No.</b>	O94892
<b>Product Type</b>	Recombinant Protein
<b>Immunogen Species</b>	Homo sapiens (Human)
<b>Purity</b>	>85% (SDS-PAGE)
<b>Sequence</b>	<p>MINAQELLTL EDVTVFTWE EWQLLGPFQK DLYRDVMLEI YSNLLSMGYQ  VSKPDALSKL ERGEEPWTME DERHSRICPE NNEVDDHLQD HLENQRMLKS  VEQYHEHNAF GNTASQTKSL CLFRENHDTF ELYIKTLKSN LSLVNQNKSC  EINNSTKFSG DGKSFLHGNY EELYSAAKFS VSTKANSTKS QVSKHQRTHE  IEKNHVCSEC GKAFVKKSQL TDHERVHTGE KPYGCTLCAK VFSRKSRLNE  HQRIHKREKS FICSECGKVF TMKSRLIEHQ RTHTGEKPYI CNECGKGFPG  KRNLIVHQRN HTGEKSYICS ECGKGF TGKS MLIHQRTHT GEKPYICSEC  GKGFTTKHYV IHHQRNHTGE KPYICNECGK GFTMKSRMIE HQRTHHTGEK  YICSECGKGF PRKSNLIVHQ RNHTVEKSYL CSECGKGF TV KSMLIHQRT  HTGEKPYTCS ECGKGFPLKS RLIVHQRTHHT GEKPYRCSEC GKGFI VNSGL  MLHQRTHHTGE KPYICNECGK GFAFKSNLVV HQRTHHTGEK FMCSECGKGF  TMKRYLIVHQ QIHTEEKSCI CSECGR GFAK ETELALHKQV HTGEKPYGCN  ECGKGF TMKS RLIVHQRTHHT GEKPFVCSEC RKAFFSSKRN LIVHQRTHNGN  KP</p>
<b>Source</b>	E.coli
<b>Target Names</b>	ZNF432
<b>Protein Names</b>	Recommended name: Zinc finger protein 432
<b>Expression Region</b>	1-652
<b>Notes</b>	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
<b>Tag Info</b>	Tag type will be determined during the manufacturing process.
<b>Protein Length</b>	Full length protein
<b>Reconstitution</b>	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
<b>Shelf Life</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.