



Recombinant Human Zinc finger protein 432 (ZNF432)

Product Code	CSB-MP026746HU
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	O94892
Product Type	Recombinant Protein
Immunogen Species	Homo sapiens (Human)
Purity	>85% (SDS-PAGE)
Sequence	<p>MINAQELLTL EDVTVEFTWE EWQLLGPFQK DLYRDVMLEI YSNLLSMGYQ VSKPDALSKL ERGEEPWTME DERHSRICPE NNEVDDHLQD HLENQRMLKS VEQYHEHNAF GNTASQTKSL CLFRENHDTF ELYIKTLKSN LSLVNQNKSC EINNSTKFSG DGKSFLHGNY EELYSAAKFS VSTKANSTKS QVSKHQRTHE IEKNHVCSEC GKAFVKKSQL TDHERVHTGE KPYGCTLCAK VFSRKSRLNE HQRIHKREKS FICSECGKVF TMKSRLIEHQ RTHTGEKPYI CNECGKGFPG KRNLIVHQRN HTGEKSYICS ECGKGFTGKS MLIHQRTHT GEKPYICSEC GKGFTTKHYV IIHQRNHTGE KPYICNECGK GFTMKSRMIE HQRTHHTGEK YICSECGKGF PRKSNLIVHQ RNHTVEKSYL CSECGKGFTV KSMLIIHQRT HTGEKPYTCS ECGKGFPLKS RLIVHQRTHHT GEKPYRCSEC GKGFI VNSGL MLHQRTHHTGE KPYICNECGK GFAFKSNLVV HQRTHHTGEK FMCSECGKGF TMKRYLIVHQ QIHTEEKSCI CSECGRFAK ETELALHKQV HTGEKPYGCN ECGKGFTMKS RLIVHQRTHHT GEKPFVCSEC RKAFFSSKRNL IVHQRTHNGN KP</p>
Source	Mammalian cell
Target Names	ZNF432
Protein Names	Recommended name: Zinc finger protein 432
Expression Region	1-652
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.