



Recombinant Human Homeobox protein SIX3 (SIX3)

Product Code	CSB-MP021350HU
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	O95343
Storage Buffer	Lyophilized from Tris/PBS-based buffer, 6% Trehalose, pH 8.0
Product Type	Recombinant Proteins
Immunogen Species	Homo sapiens (Human)
Purity	≥85% (SDS-PAGE)
Sequence	MVFRSPLDLY SSHFLLPNFA DSHHRSILLA SSGGGNGAGG GGGAGGGSGG GNGAGGGGAG GAGGGGGGGS RAPPEELSMF QLPTLNFSPE QVASVCETLE ETGDIERLGR FLWSLPVAPG ACEAINKHES ILRARAVVAF HTGNFRDLYH ILENHKFTKE SHGKLQAMWL EAHYQEA EKL RGRPLGPVDK YRVRKKFPLP RTIWDGEQKT HCFKERTRSL LREWYLQDPY PNPSKKRELA QATGLTPTQV GNWFKNRRQR DRAAAAKNRL QHQAIGPSGM RSLAEPGCPT HGSAESPSTA ASPTTSVSSL TERADTGTSI LSVTSSDSEC DV
Source	Mammalian cell
Target Names	SIX3
Protein Names	Recommended name: Homeobox protein SIX3 Alternative name(s): Sine oculis homeobox homolog 3
Expression Region	1-332
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.