



# Recombinant Human E3 ubiquitin-protein ligase ARIH2 (ARIH2)

<b>Product Code</b>	CSB-MP002069HU
<b>Storage</b>	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
<b>Uniprot No.</b>	O95376
<b>Product Type</b>	Recombinant Protein
<b>Immunogen Species</b>	Homo sapiens (Human)
<b>Purity</b>	>85% (SDS-PAGE)
<b>Sequence</b>	MSVDMNSQGS DSNEEDYDPN CEEEEEEED DPGDIEDYYV GVASDVEQQG ADAFDPEEYQ FTCLTYKESE GALNEHMTSL ASVLKVSHSV AKLILVNFHW QVSEILDYK SNSAQLLEVA RVQPNPSKHV PTSHPPHCA VCMQFVRKEN LLSLACQHQF CRSCWEQHCS VLVKDGVGVG VSCMAQDCPL RTPEDFVFPL LPNEELREKY RRYLFRDYVE SHYQLQLCPG ADCPMVIRVQ EPRARRVQCN RCNEVFCFKC RQMYHAPTDC ATIRKWLTKC ADDSETANYI SAHTKDCPKC NICIEKNGGC NHMQCSKCKH DFCWMCLGDW KTHGSEYYEC SRYKENPDIV NQSQQAQARE ALKLYLFYFE RWENHNKSLQ LEAQTYQRIH EKIQERVMNN LGTWIDWQYL QNAAKLLAKC RYTLQYTYPY AYYMESGPRK KLFYEQQAQL EAEIENLSWK VERADSYDRG DLENQMHIIE QRRRTLLKDF HDT
<b>Source</b>	Mammalian cell
<b>Target Names</b>	ARIH2
<b>Protein Names</b>	Recommended name: E3 ubiquitin-protein ligase ARIH2 Short name= ARI-2 Short name= Protein ariadne-2 homolog EC= 6.3.2.- Alternative name(s): Triad1 protein
<b>Expression Region</b>	1-493
<b>Notes</b>	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
<b>Tag Info</b>	Tag type will be determined during the manufacturing process.
<b>Protein Length</b>	Full length protein
<b>Reconstitution</b>	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
<b>Shelf Life</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.