



Recombinant Human E3 ubiquitin-protein ligase ARIH2 (ARIH2)

Product Code	CSB-YP002069HU
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	O95376
Product Type	Recombinant Protein
Immunogen Species	Homo sapiens (Human)
Purity	>85% (SDS-PAGE)
Sequence	MSVDMNSQGS DSNEEDYDPN CEEEEEEED DPGDIEDYYV GVASDVEQQG ADAFDPEEYQ FTCLTYKESE GALNEHMTSL ASVLKVSHSV AKLILVNFHW QVSEILDYK SNSAQLLEVA RVQPNPSKHV PTSHPPHCA VCMQFVRKEN LLSLACQHGF CRSCWEQHCS VLVKDGVGVG VSCMAQDCPL RTPEDFVFPL LPNEELREKY RRYLFRDYVE SHYQLQLCPG ADCPMVIRVQ EPRARRVQCN RCNEVFCFKC RQMYHAPTDC ATIRKWLTKC ADDSETANYI SAHTKDCPKC NICIEKNGGC NHMQCSKCKH DFCWMCLGDW KTHGSEYYEC SRYKENPDIV NQSQQAQARE ALKLYLFYFE RWENHNKSLQ LEAQTYQRIH EKIQERVMNN LGTWIDWQYL QNAAKLLAKC RYTLQYTYPY AYYMESGPRK KLFYEQQAQL EAEIENLSWK VERADSYDRG DLENQMHAIE QRRRTLLKDF HDT
Source	Yeast
Target Names	ARIH2
Protein Names	Recommended name: E3 ubiquitin-protein ligase ARIH2 Short name= ARI-2 Short name= Protein ariadne-2 homolog EC= 6.3.2.- Alternative name(s): Triad1 protein
Expression Region	1-493
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.