



# Recombinant Human Eppin (SPINLW1)

<b>Product Code</b>	CSB-BP022583HU
<b>Storage</b>	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
<b>Uniprot No.</b>	O95925
<b>Product Type</b>	Recombinant Protein
<b>Immunogen Species</b>	Homo sapiens (Human)
<b>Purity</b>	≥85% (SDS-PAGE)
<b>Sequence</b>	PGLTDWLFP RRCPKIREEC EFQERDVCTK DRQCQDNKKC CVFSCGKKCL DLKQDVCEMP KETGPCLAYF LHWYDYDKKDN TCSMFVYGGC QGNNNNFQSK ANCLNTCKNK RFP
<b>Source</b>	Baculovirus
<b>Target Names</b>	EPPIN
<b>Protein Names</b>	Recommended name: Eppin Alternative name(s): Cancer/testis antigen 71 Short name= CT71 Epididymal protease inhibitor Protease inhibitor WAP7 Serine protease inhibitor-like with Kunitz and WAP domains 1 WAP four- disulfide cor
<b>Expression Region</b>	22-133
<b>Notes</b>	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
<b>Tag Info</b>	Tag type will be determined during the manufacturing process.
<b>Protein Length</b>	Full Length of Mature Protein
<b>Target Details</b>	This gene encodes an epididymal protease inhibitor, which contains both kunitz-type and WAP-type four-disulfide core (WFDC) protease inhibitor consensus sequences. Most WFDC genes are localized to chromosome 20q12-q13 in two clusters: centromeric and telomeric. This gene is a member of the WFDC gene family and belongs to the telomeric cluster. Alternatively spliced transcript variants encoding distinct isoforms have been found for this gene.
<b>Reconstitution</b>	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
<b>Shelf Life</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.