



Recombinant Human TCF3 fusion partner (TFPT)

Product Code	CSB-MP023439HU
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	P0C1Z6
Product Type	Recombinant Protein
Immunogen Species	Homo sapiens (Human)
Purity	≥85% (SDS-PAGE)
Sequence	MELEQREGTM AAVGFEEFSA PPGSELALPP LFGGHILESE LETEVEFVSG GLGGSLRER DEEEEAARGR RRRQRELNRR KYQALGRRRCR EIEQVNERVL NRLHQVQRIT RRLQQERRFL MRVLDSYGDD YRASQFTIVL EDEGSQGTDA PTPGNAENEP PEKETLSPPR RTPAPPEPGS PAPGEGPSGR KRRRVPRDGR RAGNALPEL APVQIKVEED FGFEADEALD SSWVSRGPDK LLPYPTLASP ASD
Source	Mammalian cell
Target Names	TFPT
Protein Names	Recommended name: TCF3 fusion partner Alternative name(s): INO80 complex subunit F Protein FB1
Expression Region	1-253
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.