



Recombinant Human Ras-related protein Rab-40A-like (RAB40AL)

Product Code	CSB-YP019203HU
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	P0C0E4
Product Type	Recombinant Protein
Immunogen Species	Homo sapiens (Human)
Purity	>85% (SDS-PAGE)
Sequence	MSAPGSPDQA YDFLLKFLLV GDRDVGKSEI LESLQDGTAE SPYSHLGGID YKTTTILLDG QRVKLLKLDWT SGQGRFCTIF RSYSRGAQGV ILVYDIANRW SFEGMDRWIK KIEEHAPGVP KILVGNRLHL AFKRQVPREQ AQAYAERLGV TFFEVSPLCN FNIIESFTEL ARIVLLRHRL NWLGRPSKVL SLQDLCCRTI VSCTPVHLVD KLPLPIALRS HLKSFMAKG LNARMMRGLS YSLTTSSTHK RSSLCKVKIV CPPQSPPKNC TRNSCKIS
Source	Yeast
Target Names	RAB40AL
Protein Names	Recommended name: Ras-related protein Rab-40A-like Alternative name(s): Ras-like GTPase
Expression Region	1-278
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.