



Recombinant Human Histone H2B type F-M (H2BFM)

Product Code	CSB-EP010103HU-B
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	P0C1H6
Product Type	Recombinant Protein
Immunogen Species	Homo sapiens (Human)
Purity	≥85% (SDS-PAGE)
Sequence	MKCHAGSLSA NLPALKVTGW SRCVASRQPI TVGVCGPKRP LWPSSEKPFN FHTLAKLHDI GLGNDFLDMT LKAQQQKQKS HATIPKCPPL PWSPTVTLCS CYLMAAASAM AEASSETTSE EGQSIQEPKE ANSTKAQKQK RRGCRGSRRR HANRRGDSFG DSFTPYFPRV LKQVHQGLSL SQEAVSVMDS MIHDILDRIA TEAGQLAHYT KRVTITSRDI QMAVRLLLPG KMGKLAEAQG TNAALRTSLC AIWQQRK
Source	E.coli
Target Names	H2BFM
Protein Names	Recommended name: Histone H2B type F-M Alternative name(s): Histone H2B.s Short name= H2B/s
Expression Region	1-257
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.