



Recombinant Arabidopsis thaliana Transcription factor TCP12 (TCP12)

Product Code	CSB-MP366860DOA
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	A0AQW4
Product Type	Recombinant Protein
Immunogen Species	Arabidopsis thaliana (Mouse-ear cress)
Purity	>85% (SDS-PAGE)
Sequence	MFPSLDTNGY DLFDPFIPHQ TTMFPSFITH IQSPNSHHHY SSPSPFFSSD FLESFDESFL INQFLLQQQD VAANVVEPW KFCKKLELKK KNEKCVDGST SQEVQWRRTV KKRDRHSKIC TAQGPRDRRM RLSLQIARKF FDLQDMLGFD KASKTIEWLF SKSKTSIKQL KERVAASEGG GKDEHLQVDE KEKDETLKLR VSKRRTKTME SSFKTKESRE RARKRARERT MAKMKMRLFE TSETISDPHQ ETREIKITNG VQLLEKENKE QEWSNTNDVH MVEYQMDSVS IIEKFLGLTS DSSSSIFGD SEECYTSLS VRGMSTPREH NTTSIATVDE EKSPISSFSL YDYLCY
Source	Mammalian cell
Target Names	TCP12
Protein Names	Recommended name: Transcription factor TCP12 Alternative name(s): Protein BRANCHED 2
Expression Region	1-356
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.