



# Recombinant Danio rerio UPF0565 protein C2orf69 homolog (zgc:153521)

<b>Product Code</b>	CSB-MP367009DIL
<b>Storage</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
<b>Uniprot No.</b>	A0JMH2
<b>Product Type</b>	Recombinant Protein
<b>Immunogen Species</b>	Danio rerio (Zebrafish) (Brachydanio rerio)
<b>Purity</b>	>85% (SDS-PAGE)
<b>Sequence</b>	ATGQS VSGTSAAG LLRLTSVPGY DHNRVNDVLL LRPAETQRR CDNETQNLQD GNDHVVFPPG DIQNFQQEMA LQPDAPWQC WSLERVGLTL AHRFPGCHIW VIRASQMYLH KFSCYQNFVE SNLFGAPEHS ADYGAIRHLR ALLGHSMQQA GLPNPLPPLS GTSTPGPLPA GFSLTIVGFS KGCVVLNQIV YELAGARADP ELQLFLDRVS ALYWLDGGHP GGSETWVTDK CALGELASSG VAVHAHVTPY EVRDPMRAWV GREHRHFIKT LEDLGARVSH KLHFEDEPAS IDNHFRVIQE F
<b>Source</b>	Mammalian cell
<b>Target Names</b>	zgc:153521
<b>Protein Names</b>	Recommended name: UPF0565 protein C2orf69 homolog
<b>Expression Region</b>	36-341
<b>Notes</b>	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
<b>Tag Info</b>	Tag type will be determined during the manufacturing process.
<b>Protein Length</b>	Full Length of Mature Protein
<b>Reconstitution</b>	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
<b>Shelf Life</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.