



Recombinant *Gramella forsetii* Fructose-1,6-bisphosphatase class 1 1 (fbp1)

Product Code	CSB-YP367537GEQ
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	A0M0J2
Product Type	Recombinant Protein
Immunogen Species	<i>Gramella forsetii</i> (strain KT0803)
Purity	>85% (SDS-PAGE)
Sequence	MERHTTLGEF ILENQKDFPY AKGELSALLS SIRLAGKVVN QQINKAGLAE ILGKAGKENV QGESQAKLDV LANEIFVSTL RNRGEICGLA SEELEDYLVF DEEMHKNAEY IVLIDPLDGS SNIDVDITVG TIFSIYRRIS KKGTPATLED FLQPGINQVA AGYLIYGTST ILVYTTGNGV NGFTFDPGIG SFFLSHSDIK FPKKGNTYSV NEGNYVHFPQ GIKKYIKWVQ ELNEEENRPF TSRYTGSLVA DFHRNMLLGG IYLYPQGTTA PKGKLRLLYE CNPMAFLAEQ AGGKATDGSR RIMELQPQEL HERAPFIGCN KEMVEKAEFF MQ
Source	Yeast
Target Names	fbp1
Protein Names	Recommended name: Fructose-1,6-bisphosphatase class 1 1 Short name=FBPase class 1 1 EC= 3.1.3.11 Alternative name(s): D-fructose-1,6-bisphosphate 1-phosphohydrolase class 1 1
Expression Region	1-332
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.